

'Setting New Standards in Test Technology'

4. <u>PROCEDURES</u>

4.1 MATERIALS REQUIRED

- 4.1.1 Buffer: 0.1 M sodium acetate buffer (pH 4.0).
- 4.1.2 Substrate Solution: *p*-Nitrophenyl-β-Glucopyranoside (10 mM).
 Weigh 301 mg of p-Nitrophenyl-β-Glucopyranoside into a 150 ml pyrex beaker and add 95 ml of 0.1 M sodium acetate buffer (pH 4.0). Stir the solution until the substrate dissolves (about 10 min). Adjust the volume to 100 ml with 0.1 M sodium acetate buffer (pH 4.0). Add 0.02 g of sodium azide to prevent microbial contamination. Store the solution in a sealed Duran bottle. Take about 10 ml into a polypropylene tube for day-to-day use. Do not routinely use from the substrate stock. Store both containers at 4°C.
- 4.1.3 Enzyme Preparation:

With a Gilson micropipette accurately dispense 1.0 ml of well suspended enzyme solution into a 100 mL volumetric flask and adjust to volume with 0.1 M sodium acetate buffer (pH 4.0). Sequentially dilute the preparation by adding 0.5 mL of enzyme solution to 4.5 mL of buffer (i.e. 10-fold), mix the solution well, and repeat this step to get a dilution suitable for assay.

Perform these dilutions with a Gilson pipette set to 0.5 mL (for the enzyme preparation) and an Eppendorf Multipette with the 12.5 mL tip (setting 5 for 2.5 mL plus setting 4 for 2 mL) to deliver the 4.5 mL of buffer.

4.2 ASSAY PROCEDURE:

- 4.2.1 Pre-incubate substrate solution (0.2 mL) in glass test-tubes (16 x 100 mm) at 40°C for 5 min.
- 4.2.2 Pre-incubate suitably diluted enzyme solution (about 5 mL) at 40°C for 5 min.
- 4.2.3 Add 0.2 ml of enzyme solution to the bottom of tubes containing substrate solution (0.2 mL) on carefully timed intervals (about 15 sec), mix the tube contents vigorously for 5 sec., and incubate the tubes (in quadruplicate) for exactly 10 min.
- 4.2.4 At the end of the incubation time, terminate the reaction by adding 3.0 ml of 2% trisodium phosphate solution (pH 12.0) with vigorous stirring.

- 4.2.5 Prepare a Reaction Blank solution (in duplicate) by adding 3.0 mL of 2% tri-sodium phosphate solution (pH 12.0) to 0.2 mL of substrate solution, mix vigorously, and then add 0.2 mL of enzyme solution with vigorous stirring.
- 4.2.6 Measure the absorbance values (400 nm) as follows:
 - zero the spectrophotometer with distilled water.
 - Measure the absorbance of the Reaction Blank against the water and record this value (for future reference only).
 - Then zero the spectrophotometer with the Reaction Blank.
 - Measure the absorbance of all other solutions (reaction and standards) against the Reaction Blank.
- 4.2.7 Calculate activity as follows:

Units / mL =
$$\frac{\Delta E_{400}}{10} \times \frac{3.4}{0.2} \times \frac{1}{18.1} \times \frac{100}{1.0} \times \text{Dilution}$$

$$= \Delta E_{400} \ge 9.39 \ge 0.39$$

where:

$$\begin{split} \Delta E_{400} &= Absorbance \ (reaction) - Absorbance \ (blank). \\ Incubation Time &= 10 \ min. \\ Total volume \ in \ cell &= 3.4 \ ml \\ Aliquot \ assayed &= 0.2 \ ml \\ E_{mM} \ of \ p\text{-nitrophenol} \ in \ 2\% \ tri-sodium \ phosphate &= 18.1 \\ Extraction \ volume \ = 1.0 \ ml \ of \ enzyme \ in \ total \ volume \ of \ 100 \ ml \ (Original \ Extract). \\ Dilution &= further \ dilution \ of \ the \ Original \ Extract. \end{split}$$

NOTE:- For a β -glucosidase solution of 40 Units/ml, the further dilution of the "Original Enzyme Extract" is 10-fold.